

PVI ALUMINIUM PRICE LIST APRIL 2022

Please send orders to: contact@mobilitycare.net.au

39 Cambro Road, Clayton, VIC 3168

Phone: (03) 9568 8383

<https://www.mobilitycare.net.au/shop/mobility-aids/ramps/>



ALUMINIUM WHEELCHAIR RAMP PRICE LIST & SPECIFICATIONS

SINGLE FOLD RAMPS

ITEM	LENGTH	WIDTH	FOLDED SIZE	WEIGHT	Load Capacity	PRICE
SFW230	61cm (2ft)	76cm (30")	61 X 40 X 12.4cm	4.5kgs	360 kg	\$295.00
SFW330	91.5cm (3ft)	76cm (30")	95 X 40 X 12.4cm	6.8kgs	360 kg	\$395.00
SFW430	122cm (4ft)	76cm (30")	126 X 40 X 12.4cm	9.1kgs	360 kg	\$495.00
SFW530	152cm (5ft)	76cm (30")	152 X 40 X 12.4cm	10.9kgs	360 kg	\$595.00
SFW630	183cm (6ft)	76cm (30")	183 X 41 X 12.4cm	15kgs	360kg	\$725.00

MULTIFOLDING RAMPS

ITEM	LENGTH	WIDTH	FOLDED SIZE	WEIGHT	Load Capacity	PRICE
WCR530	152cm (5ft)	76cm (30")	82 X 41 X 23cm	13.6kgs	360 kg	\$745.00
WCR630	183cm (6ft)	76cm (30")	97 X 41 X 23cm	15.9kgs	360 kg	\$875.00
WCR730	213cm (7ft)	76cm (30")	112.5 X 41 X 23cm	17.7kgs	360 kg	\$995.00
WCR830	244cm (8ft)	76cm (30")	129 X 41 X 24cm	23.6kgs	360 kg	\$1,250.00

UTILITY MULTIFOLDING RAMPS

ITEM	LENGTH	WIDTH	FOLDED SIZE	WEIGHT	Load Capacity	PRICE
UTW630	183cm (6ft)	76cm (30")	97 X 41 X 27.4cm	16.3kgs	360 kg	\$955.00
UTW730	213cm (7ft)	76cm (30")	112.4 X 41 X 27.4cm	18.6kgs	360 kg	\$1,075.00
UTW830	244cm (8ft)	76cm (30")	129 X 41 X 28cm	23.6kgs	360 kg	\$1,330.00

THRESHOLD RAMPS

ITEM	LENGTH	WIDTH	Recommended Threshold Height	WEIGHT	Load Capacity	PRICE
TH1032	25.4cm (10")	81.5cm (32")	13mm - 25mm	1.8kg	136kg per axle	\$99.00
TH1232	30.4cm (12")	81.5cm (32")	19mm - 38mm	2.2kg	136kg per axle	\$125.00
TH1632	40.6cm (16")	81.5cm (32")	38mm - 51mm	2.7kg	136kg per axle	\$165.00
TH1636	40.6cm (16")	91.5cm (36")	38mm - 51mm	2.8kg	136kg per axle	\$195.00
TH2432	61cm (2ft)	81.5cm (32")	51mm - 76mm	5.4kg	136kg per axle	\$295.00
TH2436	61cm (2ft)	91.5cm (36")	51mm - 76mm	6.8kgs	136kg per axle	\$350.00

Self Supporting Threshold Ramps

ITEM	LENGTH	WIDTH	Threshold Height	Weight	Load Capacity	PRICE
SSTH1636-1.5	40.6cm (16")	91.5cm (36")	38mm or 1.5"	3.63kg	136kg per axle	\$195.00
SSTH2436-2	61cm (2ft)	91.5cm (36")	51mm or 2"	5.45kg	136kg per axle	\$250.00

ELEV8 Adjustable Self Supporting Threshold Ramps

ITEM	LENGTH	WIDTH	Threshold Height Range	Weight	Load Capacity	PRICE
ATH1232	30.4cm (12")	81.5cm (32")	25.4 - 50.8 mm (1 -2")	2.8kg	136kg per axle	\$195.00
ATH1236	30.4cm (12")	91.5cm (36")	25.4 - 50.8 mm (1 -2")	3.64 kg	136kg per axle	\$250.00
ATH2432	61cm (2ft)	81.5cm (32")	50.8 - 88.9 mm (2 - 3 1/2")	5.5 kg	136kg per axle	\$295.00
ATH2436	61cm (2ft)	91.5cm (36")	50.8 - 88.9 mm (2 - 3 1/2")	6.4 kg	136kg per axle	\$350.00
ATH3632	91.5cm (3ft)	81.5cm (32")	89 - 126.9 mm (3.5 - 5")	8.63 kg	136kg per axle	\$445.00
ATH3636	91.5cm (3ft)	91.5cm (36")	89 - 126.9 mm (3.5 - 5")	9.5 kg	136kg per axle	\$495.00
ATH4832	122cm (4ft)	81.5cm (32")	101.6 - 152.4 mm (4 - 6")	11.81 kg	136kg per axle	\$550.00
ATH4836	122cm (4ft)	91.5cm (36")	101.6 - 152.4 mm (4 - 6")	12.73 kg	136kg per axle	\$600.00

VAN RAMPS

ITEM	LENGTH	WIDTH	FOLDED HEIGHT	WEIGHT	Load Capacity	PRICE
VAN 730	213cm (7ft)	76cm (30")	91.5cm	35 kg	360 kg	\$2,195.00
VAN 830	244cm (8ft)	76cm (30")	101.6cm	39 kg	360 kg	\$2,395.00

Minimum requirement for floor to ceiling Van 730 is 94cm & Van 830 is 104.5cm

SOLID RAMPS

ITEM	LENGTH	WIDTH	WEIGHT	Load Capacity	PRICE
SL330	91.5cm (3ft)	76cm (30")	6.8kgs	360 kg	\$425.00
SL336	91.5cm (3ft)	91.5cm (36")	8.1kgs	360 kg	\$525.00
SL430	122cm (4ft)	76cm (30")	9kgs	360 kg	\$550.00
SL436	122cm (4ft)	91.5cm (36")	10.4kgs	360 kg	\$659.00
SL530	152cm (5ft)	76cm (30")	11.36kgs	360 kg	\$670.00
SL536	152cm (5ft)	91.5cm (36")	13.6kg	360 kg	\$770.00

ONTRAC SOLID RAMPS WITH HANDRAILS

ITEM	LENGTH	WIDTH	WEIGHT	Load Capacity	PRICE
XPS436	122cm (4ft)	91.5cm (36")	22.3 kg	385 kg	\$1,495.00
XPS536	152cm (5ft)	91.5cm (36")	26.82 kg	385 kg	\$1,595.00
XPS636	183cm (6ft)	91.5cm (36")	32.3 kg	385 kg	\$1,895.00
XPS736	213cm (7ft)	91.5cm (36")	36.82 kg	385 kg	\$2,095.00
XPS836	244cm (8ft)	91.5cm (36")	41.36 kg	385 kg	\$2,395.00
XPS1036	305cm (10ft)	91.5cm (36")	50 kg	385 kg	\$2,595.00

WHEEL-A-BOUT

ITEM	LENGTH	WIDTH	FOLDED SIZE	WEIGHT	WEIGHT CAPACITY	PRICE
WAB1030	305cm (10ft)	76cm (30")	168 x 94 x 41 cm	39 kgs	360 kgs	\$2,795.00
WABU1030	305cm (10ft)	76cm (30")	161 x 94 x 41 cm	45 kgs	360 kgs	\$2,795.00
WAB1230	366 cm (12ft)	76cm (30")	199 x 94 x 41 cm	38 kgs	360 kgs	\$3,095.00
WABU1230	366 cm (12ft)	76cm (30")	191 x 94 x 41 cm	45 kgs	360 kgs	\$3,095.00

Replacement Traction Tape 356 mm (14") PRICE PER 12" (FOOT) \$9.50

Prices and specifications are subject to change without notice. Maximum load capacity stated is over 2 axels, Prices are GST exclusive.