

# The Halo Bracket

**autochair**  
Innovating Vehicle Adaptations since 1983

**45,000**  
satisfied  
users  
worldwide



**Fits in two minutes**



**More compact design**



**Easy access to batteries**

## A versatile lifting bracket for small to medium scooters

The Halo Bracket has been designed for lifting small to medium sized mobility devices and provides much better access to the battery & seat post thanks to its innovative curved design. The unique shape allows scooters to be folded down up to 19% smaller than before, allowing more space in the boot or fitting into vehicles with even smaller boot areas.

**Can be fitted in just two minutes!**



Fitted in two minutes



Lifts small to medium mobility devices



19% more compact over standard brackets



Infinitely adjustable balancing point



Easy to access batteries and seat posts














Transferable from one scooter to another



Visit the Autochair website [www.autochair.co.uk](http://www.autochair.co.uk)


Autochair Ltd. Wood Street North, Meadow Lane Industrial Estate, Alfreton, Derbyshire DE55 7JR

## Typical Halo Scooter Lifting Bracket Examples\*


Technical Specifications			Maximum Lifting Weight (S.W.L.)				Compatible Hoists					
			75kg				40/80kg		80/120kg		Olympian	
Scooters the Bracket is designed for												
Range			Pride Go Go									
Scooter												
	ES 2		Ultra X		Elite Traveller		ET Plus		Sport		ET LX w/CTS	
# Wheels	3		3	4	3	4	3	4	3	4	3	4
Total Weight	32		43	45	48	50	52	53	55	58	52	53
Seat Post	Shape											
	Size		30		38		38		50		38	
Suitable Hoists	40/80	40	✓		✗		✗		✗		✗	
		80	✓		✓		✓		✓		✓	
	80/120	80	✓		✓		✓		✓		✓	
		120	✓		✓		✓		✓		✓	
	Olympian	100	✓		✓		✓		✓		✓	
		150	✓		✓		✓		✓		✓	
		200	✓		✓		✓		✓		✓	

\*The Scooter models shown have all been tested and approved, for installation of the Halo Bracket by the manufacturer.



### Simple installation









### Multiple seat post styles



### Adjustable balancing point



### Technical Data:

-  Suitable for 30mm square, 38mm and 50mm round seat posts with a minimum wall thickness of 3mm
-  Maximum Safe Working Load of 75kg
-  Suitable for use with the Mini Hoist and Olympian Hoist
-  Total unit weight of just 3kg
-  Total unit dimensions of 450x300x70mm (L,W,H)
-  Part Number 9532063

Visit the Autochair website [www.autochair.co.uk](http://www.autochair.co.uk)

Autochair Ltd. Wood Street North, Meadow Lane Industrial Estate, Alfreton, Derbyshire DE55 7JR