

Smart Transfer Person Lift

45,000
satisfied
users
worldwide

Simple
two button
operation



Get yourself mobile with The Smart Transfer

The Smart Transfer Person Lift is designed to smoothly transfer users from a wheelchair directly into a standard car seat, meaning there's no need to change your vehicle. Strong, versatile and easy to detach when not in use, the Smart Transfer Person Lift can also be used with the Smartbase for use in the home, on holiday or even in the office.



Lifts people onto the standard car seat in virtually any make or model of vehicle



Simple to install, assemble and operate



No extra lifting required



Can lift up to 23.5 stone (150kgs)



Can be used in and around the home with the Autochair Smartbase, ask for details

MPVs



4x4s



Estates



Vans



Saloons



Convertibles



How The Smart Transfer works



1 Place the sling on the seat



2 Securely attach the sling to the four point lifting attachment



3 Press the up button and then guide the person into the car



4 Unclip the sling and stow the lift ready for travel

Why thousands of people choose The Smart Transfer

- Can be used for the drivers' seat, or for front and rear seat passengers
- Simple and easy to operate
- Easily transferred to your next vehicle
- Lifting capacity of 23.5 stone (150kgs)
- Manufactured in the UK
- XS-XL slings sizes available

Dealer Details

MOBILITYCARE
mobility and rehab since 1992

For all enquiries:

PHONE: (03) 9568 8383

ADDRESS: 1/12 Marriott St, Oakleigh VIC 3166

EMAIL: info@mobilitycare.net.au

WEB: www.mobilitycare.net.au

autochair
Innovating Vehicle Adaptations since 1983

Visit the Autochair website www.autochair.co.uk

Autochair Ltd. Wood Street North, Meadow Lane Industrial Estate, Alfreton, Derbyshire DE55 7JR



Standard
installation
takes just
3-4 hours

Sling options



Standard



Double amputee



Fur lined



Mesh