

# BENOIT SYSTEMES

*Autonomy is power*



## LIGHT RANGE



[WWW.BENOITSYSTEMES.COM](http://WWW.BENOITSYSTEMES.COM)

# POWER

Using the inherent strengths of manual wheelchairs, our power add-ons give people the freedom to tackle many terrains, remove the logistical problems and turn the journeys into more about the destination than the challenge of getting there.

## THE FEATURES OF OUR POWER ADD-ONS



Tackles almost any terrain




Suitable for all wheelchairs\*



Attached in seconds



Lightweight & transportable



Discreet & unobtrusive design



Indoors / outdoors



Allowed on airplanes



Third party assistance

\*ask us for complex wheelchairs



Robust and proven technology with strong backup.



BENOIT SYSTEMES, since 1994



# HOW DOES IT WORK ?

You will be surprised by the lightness and simplicity of use of our power add-ons. They will not change either the overall dimensions of the chair, the folding ability nor the center of balance.

**1**

YOUR WHEELCHAIR + HOOKS

**2**

PLACE THE POWER ADD-ON

**3**

ADD BATTERY AND ANTI-TIP FOOT

Set-up in 15 seconds

## LIGHT DRIVE<sup>2</sup>

LIGHT DRIVE TURNS YOUR MANUAL WHEELCHAIR INTO AN ULTRA-LIGHT POWERCHAIR IN SECONDS.



DRIVE WITH A JOYSTICK!

## LIGHT ASSIST<sup>2</sup>

ATTENDANT CONTROLLED LIGHT ASSIST TAKES ALL THE EFFORT OUT OF PUSHING THE WHEELCHAIR.



PUSH THE THROTTLE HANDLE AND GO!



## THE LIGHT **DRIVE<sup>2</sup>**



**USER**  
operated

Turn your manual wheelchair into an ultra-light power chair in seconds

Transports easily in any vehicle

Free wheel when you need

Independence to go where you want

**WHEELCHAIR POWER**  
ASSISTANCE MADE EASY



## THE LIGHT **ASSIST<sup>2</sup>**



**COMPANION**  
operated

Effortless pushing of the wheelchair for the attendant

Efficient and safe maneuvers

Speed control and braking aid

Journeys together become easier

**A POWERFUL & COMPACT**  
ATTENDANT CONTROLLED POWER ADD-ON

# LIGHT DRIVE<sup>2</sup> Mini

A SMALLER VERSION FOR CHILDREN



**LIGHT DRIVE<sup>2</sup> Mini** is instinctive, simple and easy to drive.

It encourages children's independence helping with the social development.

**LIGHT DRIVE<sup>2</sup> Mini** can be refitted as the child grows allowing to change chair without changing power add-on.

# HEAVY DUTY VERSION FOR USERS ABOVE 120 KG

# LIGHT DRIVE<sup>2</sup> Plus

The extra-solid **LIGHT DRIVE<sup>2</sup> Plus** supports user's weight up to 240 kg.

Rugged and robust it enables users to enjoy the advantages of four powerful motors while keeping the lightness and mobility of a manual wheelchair.

Autonomy of at least 20 km is provided by a 18 Ah battery.



# TECHNICAL DETAILS

Our power add-ons can be fitted on almost any wheelchair starting with 32 cm seat width.

	Maximum user weight	Power add-on weight
<b>LIGHT DRIVE<sup>2</sup></b>	120 kg (265 lbs)	8,8 kg (19,4 lbs)
<b>LIGHT ASSIST<sup>2</sup></b>	120 kg (265 lbs)	8,1 kg (17,8 lbs)
<b>LIGHT DRIVE<sup>2</sup> Mini</b>	80 kg (176 lbs)	8,8 kg (19,4 lbs)
<b>LIGHT DRIVE<sup>2</sup> Plus</b>	240 kg (529 lbs)	16 kg (35,2 lbs)

- Max speed** ————— 10 km/h programmable
- Motor** ————— 2 x 130 W motors with brakes
- Motor LIGHT DRIVE<sup>2</sup> Plus** ————— 4 x 130 W motors with brakes
- Climbing ability** ————— 20 % gradient based on 70 kg person
- Manufacturer's warranty** ————— 24 months for the power add-on and the battery



OUR POWER ADD-ON IS AVAILABLE IN DIFFERENT COLORS PRESENTED ON OUR WEBSITE



12AH LITHIUM BATTERY IS COMPATIBLE WITH APPLICABLE IATA REGULATIONS FOR AIR TRAVEL

Model	Weight	Range*
— Lithium Battery	— 3,85 kg (8,5 lbs)	— 12 to 16 km (7,4 to 9,9 m)

*\*Range depends on terrain and weight of the user*



For more information visit our website [www.benoitsystemes.com](http://www.benoitsystemes.com)

# DISTRIBUTORS

## ARGENTINA

Outlandish  
www.outlandish.com.ar

## AUSTRALIA

Mobility Care  
www.mobilitycare.net.au

## BELGIUM

Mobility By Olivier  
www.mobilitybyolivier.be

## CANADA

49 Bespoke  
www.49bespoke.com

## CYPRUS

Medicom Services  
www.medicomstore.com

## DENMARK

Sunrise Medical  
www.sunrisemedical.dk

## FINLAND

Movetta  
www.movetta.fi

## IRELAND

MMS Medical  
www.mmsmedical.ie

## ITALY

Progettiamo Autonomia  
www.my-pa.it

## LUXEMBOURG

Megamed  
www.megamed.lu

## NETHERLANDS

Mobility & You  
www.mobility-you.nl

## PORTUGAL

Rehapoint  
www.rehapoint.pt

## SINGAPORE

Teng Chiang  
Holdings  
www.tchold.com

## SOUTH AFRICA

Sitwell  
www.sitwell.co.za

## SWITZERLAND

Degonda Rehab  
www.turbo-twist.ch

## SWITZERLAND

Gloor Rehabilitation  
www.gloorrehab.com

## UNITED ARAB

EMIRATES  
NMC Trading  
www.nmctrading.com

## UNITED KINGDOM

Recare  
www.recare.co.uk



WWW.BENOITSYSTEMES.COM



Follow us!

7 Rue du Pont,  
21450 Billy Les Chanceaux FRANCE

(+33) 3 80 96 51 25

export@benoitssystemes.com

CE Compliant with the applicable EU Medical  
Device Regulation 2017/745

YOUR LOCAL DEALER