



# THE LIGHT DRIVE



**POWER ADD-ON**



**THE CONTROLLER**



**AVAILABLE IN  
8 COLOURS**



**Turn your manual wheelchair into  
a powerful, ultra-light powerchair in seconds!**

**LIGHTWEIGHT, VERSATILE, QUIET, COMPACT  
& TRANSPORTS EASILY IN YOUR CAR**

# HOW DOES IT WORK?

## Wheelchair power assistance made easy with Light Drive.

- ✓ Turn your manual wheelchair into an ultra-light powerchair in seconds!
- ✓ Journeys become easier
- ✓ Tackles almost any terrain
- ✓ Transports easily in car
- ✓ Independence to go where you want



Place the power add-on onto the hooks of your wheelchair. Suitable for almost all wheelchairs from 36cm between wheels.



Slide the battery rack on the support, the connection is automatic.



The controller mounts onto a bracket which attaches to your wheelchair.



Two motors drive the wheels.



The initial balance of the wheelchair is kept, a removable anti-tip foot guarantees stability.

### Specifications

Weight (unit):	8.8kg (heaviest part)
Max weight (occupant):	120kg (depending on terrain)
Range:	12-16km
Speed:	10kph
Climbing ability:	20% gradient based on 70kg person
Battery:	Lithium ion battery (accepted for air travel by IATA)
Brakes:	Motor also acts as a downhill brake
Colour:	8 colours available

## Robust and proven technology with strong Australian backup.

***Benoit, France since 1994. MobilityCare since 1992***

YOUR LOCAL DEALER:



# MOBILITYCARE

### CONTACT

1/12 Marriott Street OAKLEIGH 3166 VIC

PHONE: (03) 9568 8383

info@mobilitycare.net.au

www.mobilitycare.net.au